



## KEY FEATURES

- Smooth, 3-Step Turret Magnification System
- LED Illumination System
- APOchromatic Lenses
- 175-400mm Working Distance Objective Lens
- Digital, Live Video and/or HD Video Capture Card System
- Widefield Eyepieces
- Straight Binocular Head (Optional: Fixed Inclined or Inclinalable)
- Green Filter
- Lifetime Warranty on Optics and Mechanics



# 935 COLPOSCOPE



FOCUS ON PRECISION.



# SEILER 935 SERIES LED COLPOSCOPE



## NEW HD VIDEO CAPTURE CARD SYSTEM



### 935 3-Step Colposcope

Small Diameter Base on Roller Casters  
with Adjustable Vertical Post Height  
3 Levels of Magnification

**Live Capture of HD Images and Video Available on All Models**



#### 0-220°

#### BINOCULAR HEAD

The inclinable head allows for a full 0-220° rotation. Seiler offers: 0-220° inclinable, straight, and 45° inclined binoculars.




#### DIGITAL DSLR CAMERA

(Sony, Canon, or Nikon)



#### ASSISTANT TEACHING HEAD

Allows a second user to view the same procedure.



#### BEAMSPLITTER

Specially designed to create the highest quality imagery. The port contains a 50/50 side for digital and 80/20 for video. Optional 50/50-50/50 available.



#### VARIO FOCUS

Never change an objective again. Stay in focus from 200mm to 350mm by simply turning the fine focus knob.



#### GO PRO HERO 5



#### MICROMANIPULATOR

The Seiler CO2 Laser Micromanipulator is intended to be used with articulated arm CO2 surgical lasers.



#### FIXED 45° HEAD FOR HRA PROCEDURES



Toll-Free: 800.489.2282 | Local: 314.968.2282 | Fax: 314.968.3601 | E-mail: [micro@seilerinst.com](mailto:micro@seilerinst.com) | [www.seilerinst.com](http://www.seilerinst.com)

MED8068 Rev.0 October 2019